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
1. AHA Guidelines 2005. *Circulation*. 2005;112 (suppl IV).
2. Mittal, S, Ayati S et al. Comparison of a rectilinear biphasic waveform with a damped sine wave monophasic waveform for transthoracic conversion of ventricular fibrillation. *JACC*, Vol 34, No 5, 1999.
3. Mittal, S, Ayati S et al. Transthoracic Cardioversion of Atrial Fibrillation Comparison of Rectilinear Biphasic Versus Damped Sine Wave Monophasic Shocks. *Circulation*. 2000;101:1282-1287.
4. W Tang, J Wang, C Young, Comparison of Rectilinear Biphasic Defibrillation To Biphasic Truncated Exponential Defibrillation For Pediatric Cardiac Life Support In A Porcine Model, Abstract presented at the Cardestim 2004.
5. Morrison I, Dorian P et al. Out-of-hospital cardiac arrest rectilinear biphasic to monophasic damped sine defibrillation waveforms with advanced life support intervention trial (ORBIT). *Resuscitation*. 2005;66:149-157.
6. van Alem et al. *Resuscitation*. 2003;58:17-24.
7. Morrison I, Dorian P et al. Out-of-hospital cardiac arrest rectilinear biphasic to monophasic damped sine defibrillation waveforms with advanced life support intervention trial (ORBIT). *Resuscitation*. 2005;66:149-157.
8. *Ibid*.
9. Niebauer MJ, Brewer JE, Chung MK and Tchou PJ. Comparison of Rectilinear Biphasic Waveform with Monophasic Damped Sine Waveform for external cardioversion of Atrial Fibrillation and Flutter. *Am J Cardiol*. 2004;93:1495-1499.
10. Li Y, Ristagno, G et al. A comparison of defibrillation efficacy between different impedance compensation techniques in high impedance porcine model. *Resuscitation*. 2009;80:1312-1317.
11. Niebauer MJ, Brewer JE et al. Rectilinear Biphasic: Defibrillation of Patients with Weight Greater Than 100 Kg. Abstract from the 7th Congress of the ERC, 2004.
12. Stiell IG, Walker RJ et al. Higher Energy Levels for Defibrillation in Out-of-Hospital Cardiac Arrest BIPHASIC Trial. A Randomized Comparison of Fixed Lower Versus Escalating. *Circulation*. 2007;115.
13. Koster RW, Walker RD, van Alem AP. Definition of successful defibrillation. *Crit. Care Med*. 2006 Dec 34 (12 Suppl).
14. Stothert JC, Hatcher TS et al. Rectilinear biphasic waveform defibrillation of out of hospital cardiac arrest *Pre-hospital Emergency Care*. 2004 Oct-Dec;8(4):388-92.

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15. Edelson DP, Abella BS, Kramer-Johansen J et al. Effects of compression depth and pre-shock pauses predict defibrillation failure during cardiac arrest. *Resuscitation*. 2006;71:137–145.

16. Kramer-Johansen, J, Myklebust H et al. Quality of out-of-hospital cardiopulmonary resuscitation with real time automated feedback: A prospective interventional study. *Resuscitation*. 2006;71:283–292.

17. White Rd, Hankins DG et al. Patient Outcomes following defibrillation with a low energy biphasic truncated exponential waveform in out of hospital cardiac arrest. *Resuscitation*, 2001;49(1):9-14.

18. Schneider T, Maartens PR et al. Multicenter, Randomized, Controlled Trial of 150J Biphasic Shocks Compared With 200- to 360J Monophasic Shocks in the Resuscitation of Out-of-Hospital Cardiac Arrest Victims. *Circulation*. 2000;102:1780-1787.

19. Koster RW, Walker RG, Chapman FW Recurrent ventricular fibrillation during advanced life support care of patients with prehospital cardiac arrest. *Resuscitation*. 2008;78:252–257.

20. Stiell IG, Walker RJ et al. Higher Energy Levels for Defibrillation in Out-of-Hospital Cardiac Arrest BIPHASIC Trial. A Randomized Comparison of Fixed Lower Versus Escalating. *Circulation*. 2007;115.

21. Koster RW, Walker RG, Chapman FW Recurrent ventricular fibrillation during advanced life support care of patients with prehospital cardiac arrest. *Resuscitation*. 2008;78:252–257.

22. Walker RG, Koster RW et al. Defibrillation probability and impedance change between shocks during resuscitation from out-of-hospital cardiac arrest. *Resuscitation*. 2009;80:773–777.

23. Koster RW, Walker RG, Chapman FW Recurrent ventricular fibrillation during advanced life support care of patients with prehospital cardiac arrest. *Resuscitation*. 2008;78:252–257.

24. Kim ML, Kim SG et al. Comparison of Rectilinear Biphasic Waveform Energy Versus Truncated Exponential Biphasic Waveform Energy for Transthoracic Cardioversion of Atrial Fibrillation *American Journal of Cardiology*, Vol. 94, December 1, 2004.

25. Alatawi F, Gurevitz O et al. Prospective, randomized comparison of two biphasic waveforms for the efficacy and safety of transthoracic biphasic cardioversion of atrial fibrillation. *Heart Rhythm*. 2005;2:382–387.

26. Neal S, Ngarmukos T et al. Comparison of the Efficacy and Safety of Two Biphasic Defibrillator Waveforms for the Conversion of Atrial Fibrillation to Sinus Rhythm. *Am J Cardiol* 2003;92:810–814.

27. Khaykin Y, Newman, D et al. Biphasic versus monophasic cardioversion in shock resistant Atrial Fibrillation. *J Cardiovasc Electrophysiology*, Vol. 14, 862-872, 2003.

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